

Product Datasheet

CD79B / CD79 Beta Antibody (Internal) (orb1537615)

Description	CD79B / CD79 Beta Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the central region of human CD79b.
Target	CD79B / CD79 Beta
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.. Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1 mg/ml
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (1:100), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.