

Product Datasheet

**C5AR2 / GPR77 / C5L2 Antibody (aa1-50)
(orb1537627)**

Description

C5AR2 / GPR77 / C5L2 Antibody (aa1-

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	A portion of amino acids 1-50 of human GPR77.
Target	C5AR2 / GPR77 / C5L2
Preservatives	PBS, 0.05% Sodium Azide. PBS, 0.05% Sodium Azide
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: LS-B1778 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B1778 was determined to be 10 ug/ml.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (0.5 - 2 µg/ml)
Expiration Date	12 months from date of receipt.