


Product Datasheet

Caspase 3 antibody (orb153764)

Description	Rabbit polyclonal antibody to Caspase 3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Caspase-3(NP_004337.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32KDa/17KDa
Uniprot ID	P42574
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200, IF/ICC: 1:20-100, IP: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...