

Product Datasheet

CD7 Antibody (clone MEM-186) (orb1537664)

Description	CD7 Antibody (clone MEM-186)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP, WB
Immunogen	Human acute myelogenous leukaemia cell line KG-1
Target	CD7
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo, IHC-P (2.5 - 5 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.