

Product Datasheet

CD44 Antibody (clone MEM-263) (orb1537677)

Description	CD44 Antibody (clone MEM-263)
Species/Host	Mouse
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, IP, WB
Immunogen	COS-7 cells (African Green Monkey)
Target	CD44
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo, IHC, IHC-P (10 µg/ml), IP, WB (2 µg/ml)
Expiration Date	12 months from date of receipt.