

Product Datasheet

Vimentin Antibody (clone VI-10) (orb1537685)

Description	Vimentin Antibody (clone VI-10)
Species/Host	Mouse
Reactivity	Gallus, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IHC-P, IP, WB
Target	Vimentin
Preservatives	TBS, pH 8, 15 mM sodium azide.. TBS, pH 8, 15 mM sodium azide.
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	ICC, IHC, IHC-P (5 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.