

Product Datasheet

MAP2 Antibody (clone MT-07) (orb1537686)

Description	MAP2 Antibody (clone MT-07)
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, IHC-Fr, IHC-P, IP, WB
Immunogen	Microtubule protein (bovine brain) enriched for kinesin
Target	MAP2
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA, ICC, IHC, IHC-Fr, IHC-P (10 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.