

Product Datasheet

SYK Antibody (aa5-360, clone SYK-01) (orb1537692)

Description	SYK Antibody (aa5-360, clone SYK-01)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IHC-P, IP, WB
Immunogen	Recombinant fragment (aa 5-360) of human Syk.
Target	SYK
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ICC, IHC, IHC-P (5 µg/ml), IP, WB (1 - 2 µg/ml)
Expiration Date	12 months from date of receipt.