

Product Datasheet

NOD2 / CARD15 Antibody (orb1537865)

Description

NOD2 / CARD15 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, IHC-P, WB
Immunogen	A synthetic peptide corresponding to 14 amino acids at the C-Terminus of human NOD2. The immunogen is located within the last 50 amino acids of NOD2.
Target	NOD2 / CARD15
Preservatives	PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C for up 3 months. Store at -20°C for up to 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: NOD2 antibody can be used for detection of Noxa by Western blot at 2 to 4 ug/mL. Antibody can also be used for immunocytochemistry starting at 10 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Western Blot: 114 kDa
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, ICC (10 µg/ml), IF, IHC, IHC-P (3 µg/ml), WB (2 - 4 µg/ml)
Expiration Date	12 months from date of receipt.