

Product Datasheet

CD95 / FAS Antibody (clone LT95) (orb1537891)

Description	CD95 / FAS Antibody (clone LT95)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P
Immunogen	HUT-78 human T cell lymphoma cell line
Target	CD95 / FAS
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo, IHC, IHC-P (10 µg/ml)
Expiration Date	12 months from date of receipt.