

Product Datasheet

ACE2 / ACE-2 Antibody (orb1538004)

Custom Reagents.

ACE2 / ACE-2 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	A synthetic peptide corresponding to amino acids near the center of human ACE2. The immunogen is located within aa 150-200 of ACE2.
Target	ACE2 / ACE-2
Preservatives	PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C for up 3 months. Store at -20°C for up to 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: ACE2 antibody can be used for the detection of ACE2 by Western blot at 0.5 to 2 ug/mL. A 90 kDa band can be detected. Antibody can also be used for immunohistochemistry starting at 2 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Western Blot:
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IF, IHC-P (2.5 µg/ml), WB (0.5 - 2 µg/ml)
Expiration Date	12 months from date of receipt.