



Product Datasheet

Collagen I antibody (orb153813)


Description	Rabbit polyclonal antibody to Collagen I
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human Collagen I/COL1A1(NP_000079.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140KDa/220KDa
Uniprot ID	P02452
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of SGC...



Immunohistochemical staining of human es...



Immunohistochemical staining of human pr...