

Product Datasheet

ITGB2 / CD18 Antibody (aa534-546, clone MEM-48, Low Endotoxin) (orb1538138)

Description	ITGB2 / CD18 Antibody (aa534-546, clone MEM-48, Low Endotoxin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC, IHC-Fr, IHC-P, IP, WB
Immunogen	Leukocytes of a patient suffering from a LGL-type leukemia
Target	ITGB2 / CD18
Preservatives	Sterile PBS, pH 7.4, Azide free. Sterile PBS, pH 7.4, Azide free
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Application notes	Does not react with: Canine
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo (2 µg/ml), Func, IHC, IHC-Fr, IHC-P (10 µg/ml), IP (2 µg/ml), WB (2 µg/ml)
Expiration Date	12 months from date of receipt.