

Product Datasheet

KRT18 / CK18 / Cytokeratin 18 Antibody (clone C-04) (orb1538139)

Description	KRT18 / CK18 / Cytokeratin 18 Antibody (clone C-04)
Species/Host	Mouse
Reactivity	Mammal
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC, IHC-P, IP, WB
Immunogen	Cytoskeleton preparation of epidermal carcinoma cell line A431
Target	KRT18 / CK18 / Cytokeratin 18
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA, Flo (1 - 2 µg/ml), ICC, IHC, IHC-P (10 µg/ml), IP (1 - 2 µg/ml), WB (1 - 2 µg/ml)
Expiration Date	12 months from date of receipt.