

Product Datasheet

HLA-E Antibody (clone MEM-E/07) (orb1538146)

Description	HLA-E Antibody (clone MEM-E/07)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, IP
Immunogen	Bacterially expressed recombinant HLA-E refolded with beta2-microglobulin and peptide
Target	HLA-E
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo (1 - 4 µg/ml), IHC, IHC-P (10 µg/ml), IP
Expiration Date	12 months from date of receipt.