

Product Datasheet

CD8A / CD8 Alpha Antibody (C-Terminus, clone P17-V) (orb1538191)

Description

CD8A / CD8 Alpha Antibody (C-

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Peptide derived from C-Terminus of human CD8.
Target	CD8A / CD8 Alpha
Preservatives	20 mM Tris-HCl, pH 8.0, 20 mg/ml rAlbumin, 0.05% Sodium Azide
Concentration	0.6 mg/ml
Storage	Store at 4°C. Do not freeze.
Note	For research use only
Application notes	Further information: Immunohistochemistry: LS-B3914 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working dilution for LS-B3914 was determined to be 1:50.
Clonality	Polyclonal
Dilution Range	IHC-Fr (1:100 - 1:200), IHC-P (40 µg/ml)
Expiration Date	12 months from date of receipt.