

Product Datasheet

KRT17 / CK17 / Cytokeratin 17 Antibody (C-Terminus, clone V21-R) (orb1538193)

Description

KRT17 / CK17 / Cytokeratin 17

	R)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-Fr, IHC-P
Immunogen	Peptide derived from C-Terminus of human cytokeratin 17.
Target	KRT17 / CK17 / Cytokeratin 17
Preservatives	20 mM Tris-HCl, pH 8.0, 20 mg/ml rAlbumin, 0.05% Sodium Azide
Concentration	0.2 mg/ml
Storage	4°C. Do not freeze.
Note	For research use only
Application notes	Further information: Immunohistochemistry: LS-B3918 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working dilution for LS-B3918 was determined to be 1:50.
Clonality	Polyclonal
Dilution Range	IHC, IHC-Fr (1:100 - 1:200), IHC-P (1:25 - 1:100)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558