

Product Datasheet

GAD65 Antibody (aa515-528) (orb1538204)

Description	GAD65 Antibody (aa515-528)
Species/Host	Goat
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Peptide with sequence C-TLEDNEERMSRLSK, from the internal region of the protein sequence according to NP_000809.1NP_001127838.1.
Target	GAD65
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (2.5 µg/ml), Peptide-ELISA (1:32000), WB (0.5 - 1.5 µg/ml)
Expiration Date	12 months from date of receipt.