

Product Datasheet

CD95 / FAS Antibody (aa224-237) (orb1538321)

Description	CD95 / FAS Antibody (aa224-237)
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Peptide with sequence KTCRKHRKENQGSH, from the internal region of the protein sequence according to NP_000034.1NP_690610.1NP_690611.1.
Target	CD95 / FAS
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Western Blot: Approx 37kD band observed in lysates of cell line MOLT4 (calculated MW of 37.7kD according to NP_000034.1). Recommended concentration: 0.1-0.3 ug/ml.
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (3.75 µg/ml), Peptide-ELISA (1:32000), WB (0.1 - 0.3 µg/ml)
Expiration Date	12 months from date of receipt.