

Product Datasheet

CLDN1 / Claudin 1 Antibody (clone 1C5-D9) (orb1538335)

Description

CLDN1 / Claudin 1 Antibody (clone

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	CLDN1 (AAH12471, 1 aa ~ 212 aa) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	CLDN1 / Claudin 1
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	1 mg/ml
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Western Blot (Cell lysate) - positive control PRC/PRF/5. Western Blot (Cell lysate) - positive control A-431. Western Blot (Transfected lysate) - positive control transfected 293T cell line. Sandwich ELISA (Recombinant protein).
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (5 µg/ml), WB
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558