

Product Datasheet

EPCAM Antibody (clone VU-1D9) (orb1538350)

Description	EPCAM Antibody (clone VU-1D9)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-Fr, IHC-P, IP, WB
Immunogen	Small cell lung carcinoma cell line H69.
Target	EPCAM
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo (1 - 3 µg/ml), ICC (4 µg/ml), IHC, IHC-Fr, IHC-P (10 µg/ml), IP, WB (1 - 4 µg/ml)
Expiration Date	12 months from date of receipt.