

Product Datasheet

IL17 Antibody (clone TC11-18H10) (orb1538444)

Description	IL17 Antibody (clone TC11-18H10)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IHC-P, WB
Immunogen	E. coli-expressed, recombinant mouse IL-17A
Target	IL17
Preservatives	BBS. BBS
Concentration	0.5 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles. Store undiluted.
Note	For research use only
Application notes	Further information: ELISA: This antibody can be used in a sandwich ELISA with LS-C41752 as a detection antibody.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA, Flo, IHC, IHC-P (10 µg/ml), WB
Expiration Date	12 months from date of receipt.