

Product Datasheet

CCDC88A / GIV / Girdin Antibody (orb1538478)

Description

CCDC88A / GIV / Girdin Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC, IHC-P, WB

Applications
Immunogen

A 16 AA synthetic peptide near the C-Terminus human Girdin. The immunogen is located within the last 50 AA of Girdin.

Target

CCDC88A / GIV / Girdin

Preservatives

PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide

Concentration

1 mg/ml

Storage

Store at 4°C for up 3 months. Store at -20°C for up to 1 year. Avoid freeze/thaw cycles.

Note

For research use only

Application notes

Further information: Girdin antibody can be used for detection of Girdin by Western blot at 1 - 2 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in rat samples and Immunofluorescence in human samples. All other applications and species not yet tested. Western Blot:

Isotype

IgG

Clonality

Polyclonal

Dilution Range

ELISA, IF, IHC, IHC-P (5 µg/ml), WB (1 - 2 µg/ml)

Expiration Date

12 months from date of receipt.