

Product Datasheet

ATP11A Antibody (orb1538532)

Description

ATP11A Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, IHC-P, WB
Immunogen	A 19 amino acid synthetic peptide near the N-Terminus of human ATP11A. The immunogen is located within the first 50 amino acids of ATP11A.
Target	ATP11A
Preservatives	PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C for up 3 months. Store at -20°C for up to 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: ATP11A antibody can be used for detection of ATP11A by Western blot at 1 ug/mL. Antibody can also be used for immunocytochemistry starting at 10 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Western Blot:
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, ICC, IF, IHC, IHC-P (5 µg/ml), WB (1 - 2 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558