

Product Datasheet

CD69 Antibody (clone H1.2F3, FITC) (orb1538543)

Description	CD69 Antibody (clone H1.2F3, FITC)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, IHC, IHC-P, IP
Immunogen	Y245 mouse dendritic epidermal T cell line
Target	CD69
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	0.5 mg/ml
Storage	Short term: 4°C; Long term: -20°C; Protect from light.
Note	For research use only
Application notes	Further information: Flow cytometry. Immunoprecipitation. In vitro T cell and NK-cell activation. Flow Cytometry: 1 ug/1E6 cells. The applications listed have been tested for the unconjugated form of this product. Other forms have not been tested.
Isotype	IgG, ₁
Clonality	Monoclonal
Dilution Range	Flo, IHC, IHC-P (10 µg/ml), IP
Expiration Date	12 months from date of receipt.