

Product Datasheet

TH / Tyrosine Hydroxylase Antibody (aa41-90) (orb1538554)

Description	TH / Tyrosine Hydroxylase Antibody (aa41-90)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P, WB
Immunogen	Synthetic peptide from human Tyrosine Hydroxylase. (AA Range: 41-90)
Target	TH / Tyrosine Hydroxylase
Preservatives	PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide, 50% Glycerol. PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide, 50% Glycerol
Concentration	1 mg/ml
Storage	Store at -20°C for up to 1 year.
Note	For research use only
Application notes	Further information: Western Blot: Predicted MW 58 kDa
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:100 - 1:500), IHC, IHC-P (1:100), Peptide-ELISA (1:20000), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.