

Product Datasheet

NEK2 Antibody (clone 2F6) (orb1538570)

Description	NEK2 Antibody (clone 2F6)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, IP, WB
Immunogen	NEK2 (AAH43502, 331 aa ~ 445 aa) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	NEK2
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Formalin-fixed paraffin-embedded sections. Sandwich ELISA: Recombinant protein. Western blot using transfected cell lysates and the recombinant protein used as the immunogen.
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (5 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.