

Product Datasheet

CD44 Antibody (clone IM7) (orb1538579)

Description	CD44 Antibody (clone IM7)
Species/Host	Rat
Reactivity	Canine, Equine, Feline, Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-Fr, IHC-P, IP, WB
Immunogen	Dexamethasone-induced cells of the SJL mouse spontaneous myeloid leukemia M1.
Target	CD44
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	Flo, ICC, IHC, IHC-Fr, IHC-P (10 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.