

Product Datasheet

TSG101 Antibody (clone 5B7) (orb1538588)

Description	TSG101 Antibody (clone 5B7)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P, WB
Immunogen	TSG101 (NP_006283, 201 aa ~ 281 aa) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	TSG101
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	1 mg/ml
Storage	Short term: -20°C; Long term: -20°C.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Formalin-fixed paraffin-embedded sections. RNAi Knockdown: Antibody validated. Sandwich ELISA: Recombinant protein. Western Blot using transfected cell lysates, cell line lysates and the recombinant protein used as the immunogen.
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	ELISA, IF (10 µg/ml), IHC, IHC-P (5 µg/ml), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.