

Product Datasheet

SIRT6 / Sirtuin 6 Antibody (aa56-105) (orb1538597)

Description	SIRT6 / Sirtuin 6 Antibody (aa56-105)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide located between aa56-105 of human SIRT6 (Q8N6T7, NP_057623). Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Monkey, Galago, Marmoset, Mouse, Hamster, Elephant, Dog, Bat, Pig, Opossum (100%); Gibbon, Rat, Bovine, Horse, Guinea pig, Platypus (92%); Turkey, Zebra finch, Chicken, Xenopus (84%).
Target	SIRT6 / Sirtuin 6
Preservatives	PBS, 0.09% sodium azide, 2% sucrose. PBS, 0.09% sodium azide, 2% sucrose
Concentration	1 mg/ml
Storage	Short term: Store at 2-8°C. Long term: Aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (1.25 µg/ml)
Expiration Date	12 months from date of receipt.