

Product Datasheet

TXNIP Antibody (C-Terminus) (orb1538650)

Description

TXNIP Antibody (C-Terminus)

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications
Immunogen

Synthetic peptide from C-Terminus of human TXNIP (Q9H3M7, NP_006463).
Percent identity by BLAST analysis:
Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Galago, Marmoset, Mouse, Rat, Sheep, Elephant, Panda, Dog, Bovine, Bat, Rabbit, Pig, Opossum, Guinea pig, Xenopus (100%).

Target

TXNIP

Preservatives

PBS, 0.09% sodium azide, 2% sucrose.
PBS, 0.09% sodium azide, 2% sucrose

Concentration

0.5 mg/ml

Storage

Short term: Store at 2-8°C. Long term:
Aliquot and store at -20°C. Avoid
freeze/thaw cycles.

Note

For research use only

Application notes

Further information: ELISA titer using peptide based assay: 1:62500. Western Blot: Suggested dilution at 1 ug/ml in 5% skim milk / PBS buffer, and HRP conjugated anti-Rabbit IgG should be diluted in 1:50000 - 100000 as second antibody.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (5 µg/ml), WB (0.2 - 1 µg/ml)

Expiration Date

12 months from date of receipt.