

Product Datasheet

CD9 Antibody (clone MM2/57) (orb1538686)

Description	CD9 Antibody (clone MM2/57)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-Fr, IHC-P, WB
Immunogen	Human platelet membrane
Target	CD9
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C or at -20°C. Store undiluted. Avoid freeze-thaw cycles. Microcentrifugation recommended if solution contains precipitate.
Note	For research use only
Application notes	Further information: Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells or 100 ul whole blood. Method sheets are available on request.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	Flo (1:100 - 1:200), IHC, IHC-Fr (1:500 - 1:1000), IHC-P (10 µg/ml), WB
Expiration Date	12 months from date of receipt.