

Product Datasheet

CD44 Antibody (clone 8E2F3) (orb1538707)

Description	CD44 Antibody (clone 8E2F3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IHC-P, WB
Immunogen	Purified recombinant fragment of human CD44 (628-699) expressed in E. coli.
Target	CD44
Preservatives	Ascitic fluid, 0.03% sodium azide. Ascitic fluid, 0.03% sodium azide
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA (1:10000), Flo (1:200 - 1:400), IF (1:200 - 1:1000), IHC, IHC-P (1:50), WB (1:100 - 1:500)
Expiration Date	12 months from date of receipt.