

Product Datasheet

CTSD / Cathepsin D Antibody (clone 3F12-1B9) (orb1538710)

Description

CTSD / Cathepsin D Antibody (clone

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, IP, WB
Immunogen	CTSD (AAH16320, 26 aa ~ 412 aa) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	CTSD / Cathepsin D
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	0.5 mg/ml
Storage	Store at -20°C or lower. Aliquot to avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Formalin-fixed paraffin-embedded sections. Sandwich ELISA: Recombinant protein. Western Blot: Transfected lysate, Recombinant protein.
Isotype	IgG1,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (5 µg/ml), IP, WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.