

Product Datasheet

CD66a / CEACAM1 Antibody (clone 2F6) (orb1538777)

Description	CD66a / CEACAM1 Antibody (clone 2F6)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant protein encoding amino acids 1-468 of human CEACAM1 isoform 6 NP_001192273.1.
Target	CD66a / CEACAM1
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	0.26 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Sandwich ELISA: Recombinant protein. Western Blot: Recombinant protein.
Isotype	IgG1,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (10 µg/ml), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.