


Product Datasheet

FLNB antibody (orb153878)

Description	Rabbit polyclonal antibody to FLNB. Filamins are a family of dimeric actin binding proteins that fun...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human FLNB(NP_001157789.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	298kDa
Uniprot ID	075369
Dilution Range	WB: 1:500-1000, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...