

Product Datasheet

ACE2 / ACE-2 Antibody (aa140-172) (orb1538835)

Description	ACE2 / ACE-2 Antibody (aa140-172)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Target	ACE2 / ACE-2
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Storage	Maintain refrigerated at 2°C to 8°C for up to 6 months. For long term storage store at -20°C.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.