

Product Datasheet

FTH1 / Ferritin Heavy Chain Antibody (clone 3F8) (orb1538912)

Description	FTH1 / Ferritin Heavy Chain Antibody (clone 3F8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant human protein purified from E. coli (HIS-Ferritin H chain).
Target	FTH1 / Ferritin Heavy Chain
Preservatives	HEPES, 0.15 M NaCl, 0.03% Sodium Azide, 50% Glycerol, 0.01% rAlbumin
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (1:200), WB (1:1000)
Expiration Date	12 months from date of receipt.