

## Product Datasheet

### KRT19 / CK19 / Cytokeratin 19 Antibody (clone RCK108) (orb1538919)

|                            |                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | KRT19 / CK19 / Cytokeratin 19 Antibody (clone RCK108)                                                       |
| <b>Species/Host</b>        | Mouse                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                |
| <b>Tested Applications</b> | FC, ICC, IHC, IHC-Fr, IHC-P, WB                                                                             |
| <b>Immunogen</b>           | Cytoskeletal preparation of human bladder carcinoma cell line T24                                           |
| <b>Target</b>              | KRT19 / CK19 / Cytokeratin 19                                                                               |
| <b>Preservatives</b>       | PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide                                                            |
| <b>Concentration</b>       | 1 mg/ml                                                                                                     |
| <b>Storage</b>             | Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles. Store undiluted. |
| <b>Note</b>                | For research use only                                                                                       |
| <b>Application notes</b>   | Further information:                                                                                        |
| <b>Isotype</b>             | IgG1                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                  |
| <b>Dilution Range</b>      | Flo (1:100 - 1:200), ICC, IHC, IHC-Fr (1:100 - 1:200), IHC-P (10 µg/ml), WB (1:100 - 1:1000)                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                             |