

Product Datasheet

HLA-DR Antibody (clone YE2/36-HLK) (orb1538930)

Description

HLA-DR Antibody (clone YE2/36-HLK)

Species/Host

Rat

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC, IHC-Fr, IHC-P

Applications
Immunogen

EHR-B Ramos cells.

Target

HLA-DR

Preservatives

PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide

Concentration

1 mg/ml

Storage

Store at 4°C or at -20°C. Store undiluted. Avoid freeze-thaw cycles. Microcentrifugation recommended if solution contains precipitate.

Note

For research use only

Application notes

Further information: Flow Cytometry: This product is routinely tested in flow cytometry on human peripheral blood lymphocytes. Use 10 ul of the suggested working dilution to label 10⁶ cells in 100ul. Immunohistochemistry: This clone has been described reacting with paraffin embedded material following PLP (periodate-lysine-paraformaldehyde) fixation (2).

Isotype

IgG2a

Clonality

Monoclonal

Dilution Range

Flo (1:50 - 1:100), IHC, IHC-Fr, IHC-P (10 µg/ml)

Expiration Date

12 months from date of receipt.