

Product Datasheet

CD79B / CD79 Beta Antibody (Intracellular, clone AT107-2) (orb1538938)

Description

CD79B / CD79 Beta Antibody

Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IHC-P, IP, WB
Immunogen	Peptide containing 20 amino acids residues from mouse CD79b conjugated to Keyhole Limpet Hemocyanin (KLH)
Target	CD79B / CD79 Beta
Preservatives	TBS, 0.09% Sodium Azide. TBS, 0.09% Sodium Azide
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Does not react with: Chicken. Further information: Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request. Membrane permeabilization is required for this application.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA, Flo (1:50 - 1:100), IHC, IHC-P (2 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.