

Product Datasheet

XAGE1E Antibody (aa58-69) (orb1538946)

Description	XAGE1E Antibody (aa58-69)
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Peptide with sequence CGFGFRRQGEDNT, from the C-Terminus of protein sequence according to NP_597673.1.
Target	XAGE1E
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (3.75 µg/ml), Peptide-ELISA (1:8000)
Expiration Date	12 months from date of receipt.