

Product Datasheet

Mitochondria Antibody (clone 113-1) (orb1539078)

Description	Mitochondria Antibody (clone 113-1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Human cell homogenate
Target	Mitochondria
Preservatives	Supplied as a liquid.. Supplied as a liquid.
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Does not react with: Mouse, Rat
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IF (1:800), IHC, IHC-P (1:50), WB (1:1000)
Expiration Date	12 months from date of receipt.