


Product Datasheet

HBXIP antibody (orb153911)

Description	Rabbit polyclonal antibody to HBXIP. HBXIP, also named as XIP, belongs to the HBXIP family. It is or...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human LAMTOR5(NP_006393.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10KDa
Uniprot ID	O43504
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...