

Product Datasheet

CTSB / Cathepsin B Antibody (aa198-227) (orb1539272)

Description	CTSB / Cathepsin B Antibody (aa198-227)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, WB
Target	CTSB / Cathepsin B
Preservatives	PBS with 0.09% (w/v) sodium azide. PBS with 0.09% (w/v) sodium azide
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2°C to 8°C for up to 6 months. For long term storage store at -20°C.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Flo (1:10 - 1:50), IHC, IHC-P (10 µg/ml), WB (1:1000)
Expiration Date	12 months from date of receipt.