



Product Datasheet

Symmetric DiMethyl-Histone H4-R3 antibody (orb153930)


Description	Rabbit polyclonal antibody to Symmetric DiMethyl-Histone H4-R3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human Symmetric DiMethyl-Histone H4-R3
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	14kDa
Uniprot ID	P62805
Entrez	8370
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...



Immunofluorescence analysis of 293T cell...



Dot-blot analysis of all sorts of methyl...