

Product Datasheet

SLC36A1 Antibody (aa1-30) (orb1539336)

Description	SLC36A1 Antibody (aa1-30)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IHC-P, WB
Target	SLC36A1
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2°C to 8°C for up to 6 months. For long term storage store at -20°C.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Flo (1:10 - 1:50), IF (1:10 - 1:50), IHC, IHC-P (5 µg/ml), WB (1:1000)
Expiration Date	12 months from date of receipt.