

Product Datasheet

A1CF / ACF Antibody (orb1539348)

Description	A1CF / ACF Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Synthetic peptide
Target	A1CF / ACF
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide. PBS, pH 7.4, 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at -20°C for up to 1 year.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA (1:20000 - 1:80000), IHC, IHC-P (5 µg/ml), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.