

Product Datasheet

PIK3CG / PI3K Gamma Antibody (orb1539353)

Description	PIK3CG / PI3K Gamma Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide derived from Human PI3K
Target	PIK3CG / PI3K Gamma
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide. PBS, pH 7.4, 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at -20°C for up to 1 year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), WB (1:2000)
Expiration Date	12 months from date of receipt.