



Product Datasheet

MonoMethyl-Histone H3-R8 antibody (orb153938)


Description	Rabbit polyclonal antibody to MonoMethyl-Histone H3-R8
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, IP, WB
Immunogen	A synthetic monomethylated peptide around R8 of human histone H3(NP_003520.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17KDa
Uniprot ID	Q16695
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:20-100, ChIP
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...



Immunofluorescence analysis of 293T cell...



Dot-blot analysis of all sorts of methyl...